

# Planmed

## DICOM Conformance Statement

### Planmed Nuance Manager Version 3.5.x

Revision 1.1  
February 26, 2013

## 1 CONFORMANCE STATEMENT OVERVIEW

Nuance Manager 3 is an application for acquiring images with Planmed Nuance or Planmed Nuance Excel mammography system. It is also able to display images received over the network or from an interchange media or a file set in local file system. Nuance Manager 3 supports sending images across the network to other systems. It supports querying a remote system for a list of DICOM objects that may then be retrieved to the local system.

**Table 1.1: NETWORK SERVICES**

<b>SOP Classes</b>	<b>User of Service (SCU)</b>	<b>Provider of Service (SCP)</b>
<b>Transfer</b>		
Digital Mammography Image Storage – For Presentation	Yes	Yes
Digital Mammography Image Storage – For Processing	Yes	Yes
Computed Radiography Image Storage	No	Yes
Grayscale Softcopy Presentation State Storage	Yes	Yes
<b>Query/Retrieve</b>		
Patient Root Query/Retrieve Information Model – FIND	Yes	No
Patient Root Query/Retrieve Information Model – MOVE	Yes	No
Study Root Query/Retrieve Information Model – FIND	Yes	No
Study Root Query/Retrieve Information Model – MOVE	Yes	No
<b>Workflow Management</b>		
Modality Worklist Information Model – FIND	Yes	No
Storage Commitment Push Model	Yes	No
Modality Performed Procedure Step	Yes	No
<b>Print Management</b>		
Basic Grayscale Print Management Meta	Yes	No
Basic Film Box	Yes	No
Basic Grayscale Image Box	Yes	No
Basic Film Session SOP Class	Yes	No
Printer	Yes	No
Print Job	Yes	No

# Planned

**Table 1.2: MEDIA SERVICES**

<b>Media Storage Application Profile</b>	<b>Write Files (FSC or FSU)</b>	<b>Read Files (FSR)</b>
<b>Compact Disc – Recordable</b>		
General Purpose CD-R	Yes	Yes
<b>DVD</b>		
General Purpose DVD-RAM	Yes	Yes

## 2 TABLE OF CONTENTS

1	CONFORMANCE STATEMENT OVERVIEW .....	2
2	TABLE OF CONTENTS.....	4
3	INTRODUCTION .....	6
3.1	REVISION HISTORY.....	6
3.2	AUDIENCE .....	6
3.3	REMARKS .....	6
3.4	TERMS AND DEFINITIONS.....	7
3.5	BASICS OF DICOM COMMUNICATION .....	9
3.6	ABBREVIATIONS .....	10
3.7	REFERENCES .....	11
4	NETWORKING .....	12
4.1	IMPLEMENTATION MODEL.....	12
4.1.1	Application data flow.....	12
4.1.2	Functional Definition of Nuance Manager 3 AE.....	13
4.1.2.1	Storage and Storage Commitment .....	13
4.1.2.2	Print.....	13
4.1.2.3	Worklist.....	13
4.1.2.4	Retrieve .....	13
4.1.2.5	MPPS .....	13
4.2	AE SPECIFICATION.....	14
4.2.1	Implementation identifying information.....	14
4.2.2	Application context name .....	14
4.2.3	Storage Specification .....	14
4.2.3.1	SOP Classes .....	14
4.2.3.2	Association Policies.....	15
4.2.3.2.1	Number of Associations.....	15
4.2.3.2.2	Asynchronous Nature.....	15
4.2.3.3	Association Initiation Policy.....	15
4.2.3.3.1	Activity – Send Image .....	15
4.2.3.3.1.1	Description and Sequencing of Activities.....	15
4.2.3.3.1.2	Proposed Presentation Contexts.....	16
4.2.3.3.1.3	SOP Specific Conformance .....	16
4.2.4	Print Management Specification.....	17
4.2.4.1	SOP Classes .....	17
4.2.4.2	Association Policies.....	17
4.2.4.2.1	General.....	17
4.2.4.2.2	Number of Associations.....	17
4.2.4.3	Association Initiation Policy.....	17
4.2.4.3.1	Print Image.....	17
4.2.4.3.1.1	Description and Sequencing of Activities.....	17
4.2.4.3.1.2	Proposed Presentation Contexts for Print image operations....	17
4.2.4.3.1.3	SOP Specific Conformance for SOP Classes .....	18

# Planned

4.2.5	Modality Worklist Specification.....	19
4.2.5.1	SOP Classes .....	19
4.2.5.2	Association Policies.....	19
4.2.5.2.1	General.....	19
4.2.5.2.2	Number of Associations.....	19
4.2.5.2.3	Asynchronous Nature.....	19
4.2.5.3	Association Initiation Policy.....	20
4.2.5.3.1	Find and Move .....	20
4.2.5.3.1.1	Description and Sequencing of Activities.....	20
4.2.5.3.1.2	Accepted Presentation Contexts .....	20
4.2.5.3.1.3	SOP Specific Conformance for SOP Class(es).....	20
4.2.6	Query/Retrieve Specification.....	21
4.2.6.1	SOP Classes .....	21
4.2.6.2	Association Policies.....	21
4.2.6.2.1	General.....	21
4.2.6.2.2	Number of Associations.....	21
4.2.6.2.3	Asynchronous Nature.....	22
4.2.6.3	Association Initiation Policy.....	22
4.2.6.3.1	Find and Move .....	22
4.2.6.3.1.1	Description and Sequencing of Activities.....	22
4.2.6.3.1.2	Proposed Presentation Contexts.....	22
4.2.6.3.1.3	SOP Specific Conformance .....	22
4.2.7	MPPS Specification .....	23
4.2.7.1	SOP Classes .....	23
4.2.7.2	Association Policies.....	23
4.2.7.2.1	General.....	23
4.2.7.2.2	Number of Associations.....	23
4.2.7.2.3	Asynchronous Nature.....	23
4.2.7.3	Association Initiation Policy.....	24
4.2.7.3.1	Image acquisition.....	24
4.2.7.3.1.1	Description and Sequencing of Activities.....	24
4.2.7.3.1.2	Proposed Presentation Contexts.....	24
4.2.7.3.1.3	SOP Specific Conformance .....	24
4.3	NETWORK INTERFACES .....	25
4.3.1	Physical Network Interface.....	25
4.3.2	Additional Protocols .....	25
4.3.3	Ipv4 and Ipv6 Support .....	25
4.4	CONFIGURATION.....	25
4.4.1	AE Title/Presentation Address Mapping .....	25
4.4.1.1	Local AE Titles.....	25
4.4.1.2	Remote AE Title/Presentation Address Mapping.....	25
5	SUPPORT OF CHARACTER SETS .....	26
6	SECURITY .....	26
7	ANNEXES.....	26
7.1	IOD CONTENTS.....	26
7.1.1	Created SOP Instance(s) .....	26

# Planned

7.1.2	Usage of Attributes form received IOD's .....	39
7.2	DATA DICTIONARY OF PRIVATE ATTRIBUTES .....	39
7.3	STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES.....	39
7.4	PRIVATE TRANSFER SYNTAXES .....	39

## 3 INTRODUCTION

### 3.1 REVISION HISTORY

Document Version	Date of Issue	Author	Reviewed by	Description
1.0	February 14, 2011	Sami Mäläskä	Erkki Lehto	First Issue
1.1	February 26, 2013	Tarek Mohsen	Erkki Lehto	Updates to Modality Worklist and Storage

### 3.2 AUDIENCE

This document is written for the people that need to understand how Nuance Manager 3 will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product. This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features.

### 3.3 REMARKS

The scope of this DICOM Conformance Statement is to facilitate integration between Nuance Manager 3 and other products supporting DICOM communication. The Conformance Statement should be read and understood in conjunction with the DICOM Standard. DICOM by itself does not guarantee interoperability. The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of different Conformance Statements is just the first step towards assessing interconnectivity and interoperability between the product and other DICOM conformant equipment.

# Planned

- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

## 3.4 TERMS AND DEFINITIONS

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

<b>Abstract Syntax</b>	the information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.
<b>Application Entity (AE)</b>	an end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.
<b>Application Entity Title</b>	the externally known name of an Application Entity, used to identify a DICOM application to other DICOM applications on the network. Application Context – the specification of the type of communication used between Application Entities. Example: DICOM network protocol.
<b>Association</b>	a network communication channel set up between Application Entities.
<b>Attribute</b>	a unit of information in an object definition; a data element identified by a tag. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).
<b>Information Object Definition (IOD)</b>	the specified set of Attributes that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The Attributes may be specified as Mandatory (Type 1), Required but possibly unknown (Type 2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.
<b>Module</b>	a set of Attributes within an Information Object Definition that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date,

# Planned

	and Patient Sex.
<b>Negotiation</b>	first phase of Association establishment that allows Application Entities to agree on the types of data to be exchanged and how that data will be encoded. Presentation Context – the set of DICOM network services used over an Association, as negotiated between Application Entities; includes Abstract Syntaxes and Transfer Syntaxes.
<b>Protocol Data Unit (PDU)</b>	a packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.
<b>Security Profile</b>	a set of mechanisms, such as encryption, user authentication, or digital signatures, used by an Application Entity to ensure confidentiality, integrity, and/or availability of exchanged DICOM data
<b>Service Class Provider (SCP)</b>	role of an Application Entity that provides a DICOM network service; typically, a server that performs operations requested by another Application Entity (Service Class User). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).
<b>Service Class User (SCU)</b>	role of an Application Entity that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)
<b>Service/Object Pair (SOP) Class</b>	the specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management. Service/Object Pair (SOP) Instance – information object; a specific occurrence of information exchanged in a SOP Class. Examples: a specific x-ray image.
<b>Tag</b>	a 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]
<b>Transfer Syntax</b>	encoding used for exchange of DICOM information objects and messages. Examples: JPEG compressed (images), little endian explicit value representation.
<b>Unique Identifier (UID)</b>	a globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object

# Planned

Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

**Value Representation (VR)** the format type of an individual DICOM data element, such as text, an integer, a person's name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.

## 3.5 BASICS OF DICOM COMMUNICATION

This section describes terminology used in this Conformance Statement for the non-specialist. The key terms used in the Conformance Statement are highlighted in *italics* below. This section is not a substitute for training about DICOM, and it makes many simplifications about the meanings of DICOM terms.

Two *Application Entities* (devices) that want to communicate with each other over a network using DICOM protocol must first agree on several things during an initial network "handshake". One of the two devices must initiate an *Association* (a connection to the other device), and ask if specific services, information, and encoding can be supported by the other device (*Negotiation*).

DICOM specifies a number of network services and types of information objects, each of which is called an *Abstract Syntax* for the Negotiation. DICOM also specifies a variety of methods for encoding data, denoted *Transfer Syntaxes*. The Negotiation allows the initiating Application Entity to propose combinations of Abstract Syntax and Transfer Syntax to be used on the Association; these combinations are called *Presentation Contexts*. The receiving Application Entity accepts the Presentation Contexts it supports.

For each Presentation Context, the Association Negotiation also allows the devices to agree on *Roles* – which one is the *Service Class User* (SCU - client) and which is the *Service Class Provider* (SCP - server). Normally the device initiating the connection is the SCU, i.e., the client system calls the server, but not always.

The Association Negotiation finally enables exchange of maximum network packet (PDU) size, security information, and network service options (called *Extended Negotiation* information).

The Application Entities, having negotiated the Association parameters, may now commence exchanging data. Common data exchanges include queries for worklists and lists of stored images, transfer of image objects and analyses (structured reports), and sending images to film printers. Each exchangeable unit of data is formatted by the sender in accordance with the appropriate *Information Object Definition*, and sent using the negotiated Transfer Syntax. There is a Default Transfer Syntax that all systems must accept, but it may not be the most efficient for some use cases. Each transfer is explicitly

# Planned

acknowledged by the receiver with a *Response Status* indicating success, failure, or that query or retrieve operations are still in process.

Two Application Entities may also communicate with each other by exchanging media (such as a CD-R). Since there is no Association Negotiation possible, they both use a *Media Application Profile* that specifies “pre-negotiated” exchange media format, Abstract Syntax, and Transfer Syntax.

## 3.6 ABBREVIATIONS

AE	Application Entity
AET	Application Entity Title
CAD	Computer Aided Detection
CD-R	Compact Disk Recordable
CR	Computed Radiography
CT	Computed Tomography
DHCP	Dynamic Host Configuration Protocol
DICOM	Digital Imaging and Communications in Medicine
DIT	Directory Information Tree (LDAP)
DN	Distinguished Name (LDAP)
DNS	Domain Name System
FSC	File-Set Creator
FSU	File-Set Updater
FSR	File-Set Reader
GSPS	Grayscale Softcopy Presentation State
HIS	Hospital Information System
HL7	Health Level 7 Standard
IHE	Integrating the Healthcare Enterprise
IOD	Information Object Definition
IPv4	Internet Protocol version 4
IPv6	Internet Protocol version 6
ISO	International Organization for Standards
LUT	Look-up Table
MG	Mammography (X-ray)
MPPS	Modality Performed Procedure Step

# Planned

MSPS	Modality Scheduled Procedure Step
MTU	Maximum Transmission Unit (IP)
MWL	Modality Worklist
O	Optional (Key Attribute)
OSI	Open Systems Interconnection
PACS	Picture Archiving and Communication System
PDU	Protocol Data Unit
R	Required (Key Attribute)
RIS	Radiology Information System.
SCP	Service Class Provider
SCU	Service Class User
SOP	Service-Object Pair
SPS	Scheduled Procedure Step
SR	Structured Reporting
TCP/IP	Transmission Control Protocol/Internet Protocol
U	Unique (Key Attribute)
UL	Upper Layer
VR	Value Representation

## **3.7 REFERENCES**


NEMA PS3	Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <a href="http://medical.nema.org/">http://medical.nema.org/</a>
----------	--

# Planned

## 4 NETWORKING

### 4.1 IMPLEMENTATION MODEL

#### 4.1.1 Application data flow



# Planmed

## 4.1.2 Functional Definition of Nuance Manager 3 AE

All communications and image transfer with remote application is accomplished utilizing the DICOM protocol over a network using the TCP/IP protocol stack.

### 4.1.2.1 Storage and Storage Commitment

Nuance Manager 3 establishes an association with a remote AE selected by the user just prior to sending a C-STORE request to that AE.

If Storage Commitment is configured to be used in Nuance Manager 3, it opens another association for sending a Storage Commitment request for the stored SOP instances to the remote AE.

Nuance Manager 3 can be configured to wait for Storage Commitment N-EVENT-REPORT in the same association that has issued the N-ACTION. Nuance Manager 3 also accepts a request for establishing an association for the N-EVENT-REPORT later.

### 4.1.2.2 Print

Nuance Manager 3 establishes an association with a remote AE selected by the user just prior to sending a Print request to that AE.

### 4.1.2.3 Worklist

Nuance Manager 3 establishes an association with a remote AE selected by the user for Modality Worklist services. When an association is requested with a SCP, Nuance Manager 3 responds with a list of SOP Class UIDs that it will accept. If a Find request is sent then it will wait for find responses.

### 4.1.2.4 Retrieve

Nuance Manager 3 establishes an association with a remote AE selected by the user for Q/R services. When an association is requested with a SCP, Nuance Manager 3 responds with a list of SOP Class UID's that it will accept. If a Find request is sent then it will wait for Find responses. If a Move request is sent, it will wait for a Move response.

### 4.1.2.5 MPPS

The Nuance Manager 3 DICOM Modality Performed Procedure Step (MPPS) SCU service is used together with DICOM Modality Worklist SCU service. If Nuance Manager 3 MPPS service is configured to be used, it will send study ID, status of study, dates, patient name in starting the exposure, and dates and complete list of images including X-ray parameters to the server after the study has been accepted.

# Planned

## 4.2 AE SPECIFICATION

### 4.2.1 Implementation identifying information

**Table 4.2-1**  
**DICOM IMPLEMENTATION CLASS AND VERSION FOR NUANCE MANAGER**

Implementation Class UID	2.16.840.1.113669.632.10.20.99.2
Implementation Version Name	NuanceManager35

### 4.2.2 Application context name

The DICOM standard application context name for DICOM 3.0 is always proposed when initiating associations:

**Table 4.2-2**  
**DICOM APPLICATION CONTEXT FOR NUANCE MANAGER**

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

### 4.2.3 Storage Specification

#### 4.2.3.1 SOP Classes

Nuance Manager 3 provides Standard Conformance to the following Storage SOP Classes:

**Table 4.2-3**  
**SOP CLASSES FOR AE STORAGE**

SOP Class Name	SOP Class UID	SCU	SCP
Digital Mammography X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.2	Yes	Yes
Digital Mammography X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	Yes	Yes
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1	No	Yes
Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.1	Yes	Yes

# Planned

## 4.2.3.2 Association Policies

### 4.2.3.2.1 Number of Associations

Table 4.2-4

NUMBER OF ASSOCIATIONS INITIATED FOR STORAGE

Maximum number of simultaneous associations	1
---	---

Table 4.2-5

NUMBER OF ASSOCIATIONS ACCEPTED FOR STORAGE

Maximum number of simultaneous associations	1
---	---

### 4.2.3.2.2 Asynchronous Nature

Nuance Manager does not support asynchronous communication (multiple outstanding transactions over a single Association).

## 4.2.3.3 Association Initiation Policy

### 4.2.3.3.1 Activity – Send Image

Nuance Manager 3 initiates a new association for the appropriate Storage Service Class that corresponds to the image requested to be transferred. The association is closed when the image has been sent to the remote DICOM network node.

#### 4.2.3.3.1.1 Description and Sequencing of Activities

A user can select studies and request them to be sent to one or more predefined destinations. Each request is forwarded to the job queue and processed individually.

When the “Auto storage” option is active, the acquired images will be forwarded to the network job queue for a pre-configured auto-send target destination(s). “Auto storage” is triggered by the Accept Study user action. Storage is invoked by the job control interface that is responsible for processing network archival tasks. The job consists of data describing the instances marked for storage and the destinations. An internal daemon process triggered by a job for a specific network destination initiates a C-STORE request to store images. If the process successfully establishes an Association to a remote Application Entity, it will transfer each marked instance one after another via the open Association. Status of the transfer is reported through the job control interface. Only one job will be active at a time. If the C-STORE Response from the remote Application contains a status other than Success or Warning, the Association is aborted. The software will retry the failed job for configured number of times (10 by default).

The Storage AE attempts to initiate a new Association in order to issue a C-STORE request. If the job contains multiple images then a new Association will be requested for each image. If the Remote AE is configured as an archive device the Storage AE will, after all images and presentation states have been sent, transmit a single Storage Commitment request (N-ACTION) in a new Association. Upon receiving the N-

# Planned

ACTION response the Storage AE will delay releasing the Association for a configurable amount of time. If no N-EVENT-REPORT is received within this time period the Association will be immediately released (i.e. notification of Storage Commitment success or failure will be received over a separate association). However, the Storage AE is capable of receiving an N-EVENT-REPORT request at any time during an association provided a Presentation Context for the Storage Commitment Push Model has been successfully negotiated (i.e. the N-ACTION is sent at the end of one association and the N-EVENT-REPORT is received during an association initiated for a subsequent send job or during an association initiated by the Remote AE for the specific purpose of sending the N-EVENT-REPORT).

#### 4.2.3.3.1.2 Proposed Presentation Contexts

The presentation contexts that are proposed by Nuance Manager 3 for the Send Image operation are specified in the following table:

**Table 4.2-6  
PRESENTATION CONTEXT TABLE**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Digital Mammography X-Ray Image Storage For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCU	None
Digital Mammography X-Ray Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.1.2	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCU	None
Computed Radiography Image Storage	1.2.840.10009.5.1.4.1.1	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCU	None
Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.1.1.1	Implicit VR Little Endian Explicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1 1.2.840.10008.1.2.2	SCU	None

All these SOP classes conform to the standard Storage Services as specified in the DICOM Standard.

#### 4.2.3.3.1.3 SOP Specific Conformance

Nuance Manager 3 sends the attributes of MG image listed in Annex 8.1. All the mandatory (type 1 and type 2) attributes are sent.

Nuance Manager 3 can also act in the role of C-STORE SCP for receiving SOP instances to be used as priors during image acquisition.

# Planned

## 4.2.4 Print Management Specification

### 4.2.4.1 SOP Classes

Nuance Manager 3 provides Standard Conformance to the following SOP Classes:

**Table 4.2-7**  
**SOP CLASSES FOR PRINT**

SOP Class Name	SOP Class UID	SCU	SCP
Basic Grayscale Print Management (META)	1.2.840.10008.5.1.1.9	Yes	No
Basic Film Session	1.2.840.10008.5.1.1.1	Yes	No
Basic Film Box	1.2.840.10008.5.1.1.2	Yes	No
Basic Grayscale Image Box	1.2.840.10008.5.1.1.4	Yes	No
Printer	1.2.840.10008.5.1.1.16	Yes	No
Print Job	1.2.840.10008.5.1.1.14	Yes	No

### 4.2.4.2 Association Policies

#### 4.2.4.2.1 General

The maximum PDU size is 65,536 bytes.

#### 4.2.4.2.2 Number of Associations

**Table 4.2-8**  
**NUMBER OF ASSOCIATIONS INITIATED FOR AE PRINT**

Maximum number of simultaneous associations	1
---	---

### 4.2.4.3 Association Initiation Policy

#### 4.2.4.3.1 Print Image

Nuance Manager 3 initiates a new association for Print Service Class. The association is closed when all images have been printed and all print jobs have completed.

##### 4.2.4.3.1.1 Description and Sequencing of Activities

Nuance Manager 3 initiates associations for the printing of images to a Basic Print SCP. Nuance Manager 3 pre-formats the images before printing.

Once the Print image association has been established, Nuance Manager 3 sends a Basic Film Session, N\_CREATE message to the Basic Print SCP. This is followed by a Basic Film Box N\_CREATE message. Nuance Manager 3 then sends a Basic Grayscale image Box, N\_SET message. Finally, an N\_ACTION message is sent to instruct the Basic Print SCP to print either as the Basic Film Session or Basic Film Box level.

##### 4.2.4.3.1.2 Proposed Presentation Contexts for Print image operations

# Planned

The presentation contexts that are proposed by Nuance Manager 3 AE for the Print Image operation are specified in the following table.

All these SOP classes conform to the standard Print Services as specified in the DICOM standard.

**Table 4.2-9  
ACCEPTABLE PRESENTATION CONTEXTS FOR PRINT**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Basic Grayscale Print Management (META)	1.2.840.10008.5.1.1.9	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Print Job	1.2.840.10008.5.1.1.14	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

### 4.2.4.3.1.3 SOP Specific Conformance for SOP Classes

Attribute values for SOP classes proposed by Nuance Manager 3 are specified in the following table.

**Table 4.2-10  
ATTRIBUTES FOR PRINT**

SOP Class Name	Command	Attribute Name	Valid Range	Default Value
Basic Film Session	N_CREATE	Number of Copies	1-10	1
		Print Priority	HIGH, MEDIUM, LOW	HIGH
		Medium Type	PAPER, CLEAR FILM, BLUE FILM	None
		Film Destination	N/A	None
		Film Session Label	N/A	None
Basic Film Session	N_ACTION	Referenced Print Job Sequence		None
Basic Film Box	N_CREATE	Image Display Format	STANDARD\1,1 STANDARD\1,2 STANDARD\2,2 STANDARD\2,3 STANDARD\3,3 STANDARD\3,4 STANDARD\3,5 STANDARD\4,4 STANDARD\4,5 STANDARD\4,6	Mandatory, no default
		Film orientation	PORTRAIT, LANDSCAPE	PORTRAIT
		Film Size ID	8INX10IN, 10INX14IN, 14INX14IN, 24CMX24CM, 10INX12IN, 11INX14IN, 14INX17IN, 24CMX30CM	None
		Magnification Type	REPLICATE, BILINEAR, CUBIC	None
		Configuration Information	This information is printer specific	None

# Planned

		Smoothing Type	N/A	None
		Border Density	BLACK,WHITE	BLACK
		Empty Image Density	BLACK,WHITE	BLACK
Basic Film Box	N_ACTION	Referenced Print Job Sequence		None
Basic Grayscale Image Box	N_SET	Image Position	1-24	Mandatory, no default
		Samples Per Pixel	1/3	None
		Photometric Interpretation	MONOCHROME1, MONOCHROME2	MONOCHROME2
		Rows	Any integer	None
		Columns	Any integer	None
		Pixel Aspect Ratio	1/1	1/1
		Pixel Representation	0 (unsigned), 1 (signed)	0
		Requested Image Size	N/A	None
Printer	N_GET/N_EVENT-REPORT	Printer Status		None
		Printer Status Info		None
		Printer Name		None
		Manufacturer		None
		Manufacturer Model Name		None
		Software Version		None

## 4.2.5 Modality Worklist Specification

### 4.2.5.1 SOP Classes

Nuance Manager 3 provides Standard Conformance to the following SOP Classes:

**Table 4.2-11**

**SOP CLASSES FOR MODALITY WORKLIST (MWL)**

SOP Class Name	SOP Class UID	SCU	SCP
Modality Worklist Find	1.2.840.10008.5.1.4.31	Yes	No

### 4.2.5.2 Association Policies

#### 4.2.5.2.1 General

The maximum PDU size is 65,536 bytes.

#### 4.2.5.2.2 Number of Associations

**Table 4.2-12**

**NUMBER OF ASSOCIATIONS INITIATED FOR MWL**

Maximum number of simultaneous associations	1
---	---

#### 4.2.5.2.3 Asynchronous Nature

Nuance Manager does not support asynchronous communication (multiple outstanding transactions over a single Association) for Modality Worklist

## 4.2.5.3 Association Initiation Policy

### 4.2.5.3.1 Find and Move

#### 4.2.5.3.1.1 Description and Sequencing of Activities

Nuance Manager 3 opens an association and performs C-FINDs. Once the association has been established, Nuance Manager 3 will send a Find message to the Modality Worklist SCP and wait for respond. The association is closed when the initiator requests that it be closed or after an error.

#### 4.2.5.3.1.2 Accepted Presentation Contexts

Worklist Management Acceptable Find execution presentation contexts for Nuance Manager 3 are:

**Table 4.2-13**  
**ACCEPTABLE PRESENTATION CONTEXTS FOR MWL**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Modality Worklist Find	1.2.840.10008.5.1.4.31	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

#### 4.2.5.3.1.3 SOP Specific Conformance for SOP Class(es)

The following attributes can be used as search criteria in C\_FIND\_RQ:

- (0008,0050) Accession Number
- (0010,0010) Patient's name
- (0010,0020) Patient ID
- (0010,0030) Patient's Birth Date
- (0040,0100) Scheduled Procedure Step Sequence
  - (0008,0060) Modality
  - (0040,0001) Scheduled Station AE Title
  - (0040,0002) Scheduled Procedure Step Start Date
- (0040,1001) Requested Procedure ID

Nuance Manager 3 reads the following attributes from a C\_FIND\_RSP message:

- (0008,0050) Accession Number
- (0008,0090) Referring Physician Name
- (0008,1110) Referenced Study Sequence
- (0010,0010) Patient's name
- (0010,0020) Patient ID
- (0010,0030) Patient's Birth Date
- (0010,0040) Patient's Sex
- (0020,000D) Study Instance UID

# Planned

- (0032,1032) Requesting Physician
- (0032,1060) Requested Procedure Description
- (0032,1064) Requested Procedure Code Sequence
- (0040,0100) Scheduled Procedure Step Sequence
  - (0008,0060) Modality
  - (0040,0001) Scheduled Station AE Title
  - (0040,0002) Scheduled Procedure Step Start Date
  - (0040,0003) Scheduled Procedure Step Start Time
  - (0040,0006) Scheduled Performing Physician's Name
  - (0040,0007) Scheduled Procedure Step Description
  - (0040,0008) Scheduled Protocol Code Sequence
  - (0040,0009) Scheduled Procedure Step ID
- (0040,1001) Requested Procedure ID
- (0040,1002) Reason For Requested Procedure
- (0040,1003) Requested Procedure Priority

## 4.2.6 Query/Retrieve Specification

### 4.2.6.1 SOP Classes

Nuance Manager 3 provides Standard Conformance to the following SOP Classes:

**Table 4.2-14**

#### **SOP CLASSES FOR QUERY/RETRIEVE**

<b>SOP Class Name</b>	<b>SOP Class UID</b>	<b>SCU</b>	<b>SCP</b>
Patient Root Query/Retrieve Information Model – Find	1.2.840.10008.5.1.4.1.2.1.1	Yes	No
Patient Root Query/Retrieve Information Model - Move	1.2.840.10008.5.1.4.1.2.1.2	Yes	No
Study Root Query/Retrieve Information Model – Find	1.2.840.10008.5.1.4.1.2.2.1	Yes	No
Study Root Query/Retrieve Information Model - Move	1.2.840.10008.5.1.4.1.2.2.2	Yes	No

### 4.2.6.2 Association Policies

#### 4.2.6.2.1 General

The maximum PDU size is 65,536 bytes.

#### 4.2.6.2.2 Number of Associations

**Table 4.2-15**

#### **NUMBER OF ASSOCIATIONS AS AN ASSOCIATION INITIATOR FOR QUERY/RETRIEVE**

Maximum number of simultaneous associations	1
---	---

# Planned

## 4.2.6.2.3 Asynchronous Nature

Nuance Manager does not support asynchronous communication (multiple outstanding transactions over a single Association) for Query/Retrieve

## 4.2.6.3 Association Initiation Policy

### 4.2.6.3.1 Find and Move

#### 4.2.6.3.1.1 Description and Sequencing of Activities

Nuance Manager 3 initiates an association for the appropriate Query/Retrieve Service Class that corresponds to the set of images requested to be transferred. Once the association has been established, Nuance Manager 3 sends Find Q/R message (C-FIND). After response has been received, Nuance Manager 3 sends a request for a Move Service (C-MOVE) and waits for an incoming Storage association. The association is closed when all queries or moves have been sent to the remote DICOM network node.

#### 4.2.6.3.1.2 Proposed Presentation Contexts

Table 4.2-16

### PROPOSED PRESENTATION CONTEXTS FOR QUERY/RETRIEVE

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Patient Root Query/Retrieve Information Model – Find	1.2.840.10008.5.1.4.1.2.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
Patient Root Query/Retrieve Information Model - Move	1.2.840.10008.5.1.4.1.2.1.2	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

Table 4.2-17

### EXTENDED NEGOTIATION AS A SCU

SOP Class Name	SOP Class UID	Extended Negotiation

#### 4.2.6.3.1.3 SOP Specific Conformance

The following attributes can be used as search criteria in C\_FIND\_RQ.

- (0010,0010) Patient's name
- (0010,0020) Patient ID
- (0010,0030) Patient's Birth Date
- (0008,0020) Study Date

Nuance Manager reads the following attributes from a C\_FIND\_RSP message:

# Planned

- (0008,0020) Study Date
- (0008,0030) Study Time
- (0008,0050) Accession Number
- (0008,0050) Accession Number
- (0008,0060) Modality
- (0008,0061) Modalities In Study
- (0008,0090) Referring Physician Name
- (0008,1030) Study Description
- (0010,0010) Patient's name
- (0010,0020) Patient ID
- (0010,0030) Patient's Birth Date
- (0020,000D) Study Instance UID
- (0020,0010) Study ID
- (0020,000E) Series Instance UID
- (0020,0011) Series Number

## 4.2.7 MPPS Specification

### 4.2.7.1 SOP Classes

Nuance Manager 3 provides Standard Conformance to the following SOP Classes:

Table 4.2-18  
SOP CLASS(ES) FOR MPPS

SOP Class Name	SOP Class UID	SCU	SCP
Modality Performed Procedure Step	1.2.840.10008.5.1.4.1.2.1.1	Yes	No

### 4.2.7.2 Association Policies

#### 4.2.7.2.1 General

The maximum PDU size is 65,536 bytes.

#### 4.2.7.2.2 Number of Associations

Table 4.2-19  
NUMBER OF ASSOCIATIONS AS AN ASSOCIATION INITIATOR FOR MPPS

Maximum number of simultaneous associations	1
---	---

#### 4.2.7.2.3 Asynchronous Nature

Nuance Manager does not support asynchronous communication (multiple outstanding transactions over a single Association) for MPPS.

# Planned

## 4.2.7.3 Association Initiation Policy

### 4.2.7.3.1 Image acquisition

#### 4.2.7.3.1.1 Description and Sequencing of Activities

Nuance Manager 3 will initiate association as a MPPS when the local operator requests to start a new study for acquire a set of images for a patient selected via Worklist.

When the association has been established Nuance Manager 3 invokes either an N-CREATE or N-SET request to the server. When starting a new study Nuance Manager 3 sends N-CREATE request to the server. When status of the MPPS instance is to be updated, Nuance Manager 3 will initiate the MPPS N-SET service request to update the status of the MPPS instance. The COMPLETE status will be finally delivered with the MPPS N-SET request after all associated images have been acquired.

#### 4.2.7.3.1.2 Proposed Presentation Contexts

Table 4.2-20

### PROPOSED PRESENTATION CONTEXTS FOR MPPS

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Modality Performed Procedure Step	1.2.840.10008.3.1.2.3.3	Implicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.2	SCU	None

#### 4.2.7.3.1.3 SOP Specific Conformance

The following attributes are provided

- (0008,0050) Accession Number
- (0008,0060) Modality
- (0008,1050) Performing Physician's Name
  
- (0010,0010) Patient's Name
- (0010,0020) Patient ID
- (0010,0030) Patient's Birth Date
- (0010,0040) Patient's Sex
  
- (0018,0060) KVp
- (0018,8151) X-Ray Tube Current in uA
- (0018,1150) Exposure Time (msec)
- (0018,1030) Protocol Name
  
- (0020,000D) Study Instance UID
- (0020,0010) Study ID
- (0020,000E) Series Instance UID

# Planned

- (0040,1001) Requested Procedure ID
- (0040,0253) Performed Procedure Step ID
- (0040,0241) Performed Station AE Title
- (0040,0242) Performed Station Name
- (0040,0244) Performed Procedure Step Start Date
- (0040,0245) Performed Procedure Step Start Time
- (0040,0252) Performed Procedure Step Status
- (0040,0250) Performed Procedure Step End Date
- (0040,0251) Performed Procedure Step End Time
- (0040,0340) Performed Series Sequence
- (0040,0301) Total Number of Exposures
- (0040,0270) Scheduled Step Attribute Sequence

## 4.3 NETWORK INTERFACES

### 4.3.1 Physical Network Interface

Nuance Manager 3 runs on Windows XP and Windows 7 platforms and utilizes their TCP/IP support. Hence it is able to use any (TCP/IP) Physical Network Interface that Windows supports.

### 4.3.2 Additional Protocols

Nuance Manager 3 supports TCP/IP protocol only.

### 4.3.3 Ipv4 and Ipv6 Support

IPv4 is supported.

IPv6 has not been tested and therefore not yet supported.

## 4.4 CONFIGURATION

### 4.4.1 AE Title/Presentation Address Mapping

Presentation address mapping is configured in Nuance Manager 3/DICOM Settings. Please see Nuance Manager Installation Manual for details.

#### 4.4.1.1 Local AE Titles

Table 4.4-1  
AE TITLE CONFIGURATION TABLE

Application Entity	Default AE Title	Default TCP/IP Port
Nuance Manager	PLANMED_AWS_1	10410

#### 4.4.1.2 Remote AE Title/Presentation Address Mapping

Please see Nuance Manager Installation Manual for details.

# Planned

## 5 SUPPORT OF CHARACTER SETS

Character set "ISO\_IR 100" is used by default.

## 6 SECURITY

Nuance Manager 3 does not implement any of the "Secure Use Profiles" defined in PS 3.15 (section 6.1 and Annex A), nor does it implement any of the "Secure Transport Connection Profiles" as defined in PS 3.15 (section 6.2 and Annex B).

Nuance Manager 3 does not implement the "Digital Signature Profile" as defined in PS 3.15 (section 6.3 and Annex C), and the "Media Storage Security Profiles" as defined in PS 3.15 (section 6.4 and Annex D) are not applicable to Nuance Manager 3.

## 7 ANNEXES

### 7.1 IOD CONTENTS

#### 7.1.1 Created SOP Instance(s)

Tables 7.1-x specify the attributes of a Digital Mammography X-Ray image transmitted by Nuance Manager 3 application.

The abbreviations used in "Presence of Module" column are:

VNAP	Value Not Always Present (attribute sent zero length if no value is present)
ANAP	Attribute Not Always Present
ALWAYS	Always Present
EMPTY	Attribute is sent without a value

The abbreviations used in the "Source" column are:

MWL	the attribute value source Modality Worklist
USER	the attribute value source is from User input
AUTO	the attribute value is generated automatically
MPPS	the attribute value is the same as that use for Modality Performed Procedure Step
CONFIG	the attribute value source is a configurable parameter

**Table 7.1-1a**  
**IOD OF CREATED MG SOP INSTANCES**

IE	Module	Reference	Presence of Module
Patient	Patient	Table 7.1-2	ALWAYS
Study	General Study	Table 7.1-3	ALWAYS

# Planned

	Patient Study	Table 7.1-4	ANAP
Series	General Series	Table 7.1-5	ALWAYS
	DX Series	Table 7.1-6	ALWAYS
	Mammography Series	Table 7.1-7	ALWAYS
	Frame of Reference	Table 7.1-8	ALWAYS
Equipment	General Equipment	Table 7.1-9	ALWAYS
Image	General Image	Table 7.1-10	ALWAYS
	Image Pixel	Table 7.1-11	ALWAYS
	DX Anatomy Imaged	Table 7.1-12	ALWAYS
	DX Image	Table 7.1-13	ALWAYS
	DX Detector	Table 7.1-14	ALWAYS
	X-Ray Acquisition Dose	Table 7.1-15	ANAP
	X-Ray Generation	Table 7.1-16	ANAP
	X-Ray Filtration	Table 7.1-17	ANAP
	X-Ray Grid	Table 7.1-18	ANAP
	Mammography Image	Table 7.1-19	ALWAYS
	VOI LUT	Table 7.1-20	VNAP
	Acquisition Context	Table 7.1-21	ALWAYS
	SOP Common	Table 7.1-22	ALWAYS

Table 7.1-1b

## IOD OF CREATED GRAYSCALE SOFTCOPY PRESENTATION STATE SOP INSTANCES

IE	Module	Reference	Presence of Module
Patient	Patient	Table 7.1-2	ALWAYS
Study	General Study	Table 7.1-3	ALWAYS
	Patient Study	Table 7.1-4	ANAP
Series	General Series	Table 7.1-5	ALWAYS
	Presentation Series	Table 7.1-23	ALWAYS
Equipment	General Equipment	Table 7.1-9	ALWAYS
Presentation State	Presentation State Identification	Table 7.1-24	ALWAYS
	Presentation State Relationship	Table 7.1-25	ALWAYS
	Displayed Area	Table 7.1-26	ALWAYS
	Graphic Annotation	Table 7.1-27	ANAP
	Graphic Layer	Table 7.1-28	ANAP
	Modality LUT	Table 7.1-29	ANAP
	Softcopy VOI LUT	Table 7.1-30	ANAP
	Softcopy Presentation LUT	Table 7.1-31	ALWAYS
SOP Common	Table 7.1-22	ALWAYS	

# Planned

**Table 7.1-2  
PATIENT MODULE OF CREATED SOP INSTANCES**

<b>Attribute Name</b>	<b>Tag</b>	<b>VR</b>	<b>Value</b>	<b>Presence of Value</b>	<b>Source</b>
Patient Name	(0010,0010)	PN	From Modality Worklist or user input. Values supplied via Modality Worklist will be entered as received. Values supplied via user input will contain all the components that the user entered. Maximum 64 characters.	ALWAYS	MWL / USER
Patient ID	(0010,0020)	LO	From Modality Worklist or user input. Maximum 64 characters.	ALWAYS	MWL / USER
Patient's Birth Date	(0010,0030)	DA	From Modality Worklist or user input	VNAP	MWL / USER
Patient's Sex	(0010,0040)	CS	From Modality Worklist or user input	VNAP	MWL / USER
Patient Comments	(0010,4000)	LT	From User Input. Maximum 1024 characters.	VNAP	USER

**Table 7.1-3  
GENERAL STUDY MODULE OF CREATED SOP INSTANCES**

<b>Attribute Name</b>	<b>Tag</b>	<b>VR</b>	<b>Value</b>	<b>Presence of Value</b>	<b>Source</b>
Study Instance UID	(0020,000D)	UI	From Modality Worklist or generated	ALWAYS	MWL or AUTO
Study Date	(0008,0020)	DA		ALWAYS	AUTO
Study Time	(0008,0030)	TM		ALWAYS	AUTO
Referring Physician's Name	(0008,0090)	PN	From Modality Worklist or user input	VNAP	MWL / USER
Study ID	(0020,0010)	SH	From Modality Worklist or generated	ALWAYS	AUTO
Accession Number	(0008,0050)	SH	Accession number for the study	VNAP	AUTO
Study Description	(0008,1030)	LO	From user input	ANAP	MWL/USER

**Table 7.1-4  
PATIENT STUDY MODULE OF CREATED SOP INSTANCES**

<b>Attribute Name</b>	<b>Tag</b>	<b>VR</b>	<b>Value</b>	<b>Presence of Value</b>	<b>Source</b>
Patient's Age	(0010,1010)	AS	From Modality Worklist or user input	ANAP	MWL/USER
Patient's Size	(0010,1020)	DS	From Modality Worklist or user input	ANAP	MWL
Patient's Weight	(0010,1030)	DS	From Modality Worklist or user input	ANAP	MWL
Occupation	(0010,2180)	SH		ANAP	MWL

**Table 7.1-5  
GENERAL SERIES MODULE OF CREATED SOP INSTANCES**

<b>Attribute Name</b>	<b>Tag</b>	<b>VR</b>	<b>Value</b>	<b>Presence of Value</b>	<b>Source</b>
Modality	(0008,0060)	CS	“PR” or “MG”	ALWAYS	AUTO
Series Instance UID	(0020,000E)	UI	Generated	ALWAYS	AUTO
Series Number	(0020,0011)	IS	Generated	ALWAYS	AUTO
Laterality	(0020,0060)	CS	If a user has chosen one image per series then this attribute is always present otherwise never.	ANAP	USER / CONFIG
Series Date	(0008,0021)	DA		ALWAYS	AUTO
Series Time	(0008,0031)	TM		ALWAYS	AUTO
Protocol Name	(0018,1030)	LO		ANAP	USER
Series Description	(0008,103E)	LO	Laterality – Projection If configuration option “Create new series for each image” is selected	ANAP	USER or AUTO
Operators' Name	(0008,1070)	PN		ANAP	MWL/USER
Referenced Performed Procedure Step Sequence	(0008,1111)	SQ	Identifies the MPPS SOP Instance to which this image is related if MPPS is used	ANAP	AUTO
>Referenced SOP Class UID	(0008,1150)	UI	MPPS SOP Class UID	ANAP	AUTO
>Referenced SOP Instance UID	(0008,1155)	UI	MPPS SOP Instance UID	ANAP	AUTO
Body Part Examined	(0018,0015)	CS	“BREAST”	ALWAYS	AUTO

**Table 7.1-6  
DX SERIES MODULE OF CREATED SOP INSTANCES**

<b>Attribute Name</b>	<b>Tag</b>	<b>VR</b>	<b>Value</b>	<b>Presence of Value</b>	<b>Source</b>
Modality	(0008,0060)	CS	“MG”	ALWAYS	AUTO
>Referenced SOP Class UID	(0008,1150)	UI	1.2.840.10008.3.1.2.3.3	ANAP	AUTO
>Referenced SOP Instance UID	(0008,1155)	UI		ANAP	AUTO
Presentation Intent Type	(0008,0068)	CS	FOR PROCESSING or FOR PRESENTATION	ALWAYS	AUTO

**Table 7.1-7  
MAMMOGRAPHY SERIES MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Modality	(0008,0060)	CS	"MG"	ALWAYS	AUTO

**Table 7.1-8  
FRAME OF REFERENCE MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Position Reference Indicator	(0020,1040)	LO	-	EMPTY	AUTO

**Table 7.1-9  
GENERAL EQUIPMENT MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Manufacturer	(0008,0070)	LO	Planmed	VNAP	AUTO
Institution Name	(0008,0080)	LO		VNAP	CONFIG
Institution Address	(0008,0081)	ST		VNAP	CONFIG
Station Name	(0008,1010)	SH		ANAP	CONFIG
Institutional Department Name	(0008,1040)	LO		ANAP	CONFIG
Manufacturer's Model Name	(0008,1090)	LO		ALWAYS	AUTO
Device Serial Number	(0018,1000)	LO		ANAP	AUTO
Software Versions	(0018,1020)	LO		ANAP	AUTO
Date of Last Calibration	(0018,1200)	DA		ANAP	AUTO
Time of Last Calibration	(0018,1201)	TM		ANAP	AUTO
Pixel Padding Value	(0028,0120)	US		ANAP	AUTO

**Table 7.1-10  
GENERAL IMAGE MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Instance Number	(0020,0013)	IS		ALWAYS	AUTO
Patient Orientation	(0020,0020)	CS		VNAP	AUTO
Content Date	(0008,0023)	DA		ALWAYS	AUTO
Content Time	(0008,0033)	TM		ALWAYS	AUTO
Image Type	(0008,0008)	CS	ORIGINAL, PRIMARY if	ANAP	AUTO

# Planned

			image is not for presentation		
Acquisition Date	(0008,0022)	DA		ALWAYS	AUTO
Acquisition Time	(0008,0032)	TM		ALWAYS	AUTO
Source Image Sequence	(0008,2112)	SQ	If FOR PRESENTATION image. This sequence refers to the original image if present.	ANAP	AUTO
> Referenced SOP Class UID	(0008,1150)	UI		ANAP	AUTO
> Referenced SOP Instance UID	(0008,1155)	UI		ANAP	AUTO
>Spatial Locations Preserved	(0028,135A)	CS		ANAP	
Image Comments	(0020,4000)	LT		ANAP	USER
Burned In Annotation	(0028,0301)	CS	"NO"	ALWAYS	USER
Lossy Image Compression	(0028,2110)	CS	"00"	ALWAYS	AUTO
Presentation LUT Shape	(2050,0020)	CS		ALWAYS	AUTO

**Table 7.1-11  
IMAGE PIXEL MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Samples per Pixel	(0028,0002)	US	1	ALWAYS	AUTO
Photometric Interpretation	(0028,0004)	CS	MONOCHROME1 or MONOCHROME2	ALWAYS	AUTO
Rows	(0028,0010)	US		ALWAYS	AUTO
Columns	(0028,0011)	US		ALWAYS	AUTO
Bits Allocated	(0028,0100)	US	16	ALWAYS	AUTO
Bits Stored	(0028,0101)	US	12 or 16	ALWAYS	AUTO
High Bit	(0028,0102)	US	11 or 15	ALWAYS	AUTO
Pixel Representation	(0028,0103)	US	0 (Unsigned)	ALWAYS	AUTO
Pixel Data	(7FE0,0010)			ALWAYS	AUTO

**Table 7.1-12  
DX ANATOMY IMAGED MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Image Laterality	(0020,0062)	CS	R or L	ALWAYS	USER at device
Anatomic Region	(0008,2218)	SQ	I	ANAP	AUTO

# Planned

Sequence					
>Code Value	(0008,0100)	SH	T-04000	ALWAYS	AUTO
>Coding Scheme Designator	(0008,0102)	SH	SNM3	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	BREAST	ALWAYS	AUTO

**Table 7.1-13  
DX IMAGE MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Image Type	(0008,0008)	CS		ALWAYS	AUTO
Samples per Pixel	(0028,0002)	US	1	ALWAYS	AUTO
Photometric Interpretation	(0028,0004)	CS	MONOCHROME1 or MONOCHROME2	ALWAYS	AUTO
Bits Allocated	(0028,0100)	US	12 or 16	ALWAYS	AUTO
Bits Stored	(0028,0101)	US	12 or 16	ALWAYS	AUTO
High Bit	(0028,0102)	US	11 or 15	ALWAYS	AUTO
Pixel Representation	(0028,0103)	US	0 (unsigned)	ALWAYS	AUTO
Pixel Intensity Relationship	(0028,1040)	CS	LIN	ALWAYS	AUTO
Pixel Intensity Relationship Sign	(0028,1041)	SS	-1	ALWAYS	AUTO
Rescale Intercept	(0028,1052)	DS	0	ALWAYS	AUTO
Rescale Slope	(0028,1053)	DS	1	ALWAYS	AUTO
Rescale Type	(0028,1054)	LO	US	ALWAYS	AUTO
Presentation LUT Shape	(2050,0020)	CS	IDENTITY or INVERSE	ALWAYS	AUTO
Lossy Image Compression	(0028,2110)	CS	00	ALWAYS	AUTO
Acquisition Device Processing Description	(0018,1400)	LO	Image processing parameter file name for FOR PRESENTATION images	VNAP	AUTO
Patient Orientation	(0020,0020)	CS		ALWAYS	AUTO
Burned In Annotation	(0028,0301)	CS	YES or NO	ALWAYS	AUTO
Window Center	(0028,1050)	DS		ALWAYS	AUTO
Window Width	(0028,1051)	DS		ALWAYS	AUTO

**Table 7.1-14  
DX DETECTOR MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Detector Type	(0018,7004)	CS		VNAP	AUTO
Detector ID	(0018,700A)	SH		VNAP	AUTO
Date of Last	(0018,700C)	DA		VNAP	AUTO

# Planned

Detector Calibration					
Time of Last Detector Calibration	(0018,700E)	TM		VNAP	AUTO
Detector Temperature	(0018,7001)	DS		VNAP	AUTO
Imager Pixel Spacing	(0018,1164)	DS		VNAP	AUTO

**Table 7.1-15  
X-RAY ACQUISITION DOSE MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
KVP	(0018,0060)	DS	From x-ray unit	ALWAYS	AUTO
X-Ray Tube Current	(0018,1151)	IS	From x-ray unit	ALWAYS	AUTO
Exposure Time	(0018,1150)	IS	From x-ray unit	ALWAYS	AUTO
Exposure	(0018,1152)	IS	From x-ray unit	ALWAYS	AUTO
Exposure in $\mu$ As	(0018,1153)	IS	From x-ray unit	ALWAYS	AUTO
Body Part Thickness	(0018,11A0)	DS	From x-ray unit	VNAP	AUTO
Relative X-Ray Exposure	(0018,1405)	IS	Entrance dose in $10 * \mu$ Gy	ALWAYS	AUTO
Entrance Dose in mGy	(0040,8302)	DS	From x-ray unit	VNAP	AUTO
Organ Dose	(0040,0316)	DS	From x-ray unit	VNAP	AUTO
Organ Exposed	(0040,0318)	CS	"BREAST"	ALWAYS	AUTO
Anode Target Material	(0018,1191)	CS	From x-ray unit	ALWAYS	AUTO
Filter Type	(0018,1160)	SH	From x-ray unit	ALWAYS	AUTO
Filter Material	(0018,7050)	CS	From x-ray unit	ALWAYS	AUTO
Grid	(0018,1166)	CS	From x-ray unit	ALWAYS	AUTO
Rectification Type	(0018,1156)	CS	"CONST POTENTIAL"	ALWAYS	AUTO

**Table 7.1-16  
X-RAY GENERATION MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
KVP	(0018,0060)	DS	From x-ray unit	ALWAYS	AUTO
X-Ray Tube Current	(0018,1151)	IS	From x-ray unit	ALWAYS	AUTO
Exposure Time	(0018,1150)	IS	From x-ray unit	ALWAYS	AUTO
Exposure	(0018,1152)	IS	From x-ray unit	ALWAYS	AUTO
Exposure in $\mu$ As	(0018,1153)	IS	From x-ray unit	ALWAYS	AUTO
Exposure Control Mode	(0018,7060)	CS	From x-ray unit	ALWAYS	AUTO
Exposure	(0018,7062)	LT	From x-ray unit	ALWAYS	AUTO

# Planned

Control Mode Description					
Focal Spot	(0018,1190)	DS	From x-ray unit	ALWAYS	AUTO
Anode Target Material	(0018,1191)	CS	From x-ray unit	ALWAYS	AUTO
Rectification Type	(0018,1156)	CS	From x-ray unit	ALWAYS	AUTO

**Table 7.1-17  
X-RAY FILTRATION MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Filter Type	(0018,1160)	SH	From x-ray unit	ALWAYS	AUTO
Filter Material	(0018,7050)	CS	From x-ray unit	ALWAYS	AUTO

**Table 7.1-18  
X-RAY GRID MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Grid	(0018,1166)	CS	From x-ray unit	ALWAYS	AUTO

**Table 7.1-19  
MAMMOGRAPHY IMAGE MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Image Type	(0008,0008)	CS	ORIGINAL or PRIMARY	ALWAYS	AUTO
Positioner Type	(0018,1508)	CS	MAMMOGRAPHY	ALWAYS	AUTO
Positioner Primary Angle	(0018,1510)	DS	From x-ray unit	ALWAYS	AUTO
Positioner Secondary Angle	(0018,1511)	DS	From x-ray unit	ANAP	AUTO
Image Laterality	(0020,0062)	CS	From x-ray unit	ALWAYS	AUTO
Organ Exposed	(0040,0318)	CS	"BREAST"	ALWAYS	AUTO
Breast Implant Present	(0028,1300)	CS	From user input	ANAP	USER
Partial View	(0028,1350)	CS	YES or NO	ALWAYS	USER
Partial View Code Sequence	(0028,1352)	SQ	If Partial View is YES	ANAP	AUTO
>Code Value	(0008,0100)	SH		VNAP	AUTO
>Code Meaning	(0008,0104)	LO		VNAP	AUTO
Anatomic Region Sequence	(0008,2218)	SQ		ALWAYS	AUTO
>Code Value	(0008,0100)	SH	T-04000	ALWAYS	AUTO
>Coding Scheme	(0008,0102)	SH	SNM3	ALWAYS	AUTO

# Planned

Designator					
>Code Meaning	(0008,0104)	LO	BREAST	ALWAYS	AUTO

**Table 7.1-20  
VOI LUT MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Window Center	(0028,1050)	DS		ALWAYS	AUTO
Window Width	(0028,1051)	DS		ALWAYS	AUTO
Window Center and Width Explanation	(0028,1055)	LT		ANAP	AUTO

**Table 7.1-21  
ACQUISITION CONTEXT MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Acquisition Context Sequence	(0040,0555)	SQ		EMPTY	AUTO

**Table 7.1-22  
SOP COMMON MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
SOP Class UID	(0008,0016)	UI	1.2.840.10008.5.1.4.1.1.1.2 or 1.2.840.10008.5.1.4.1.1.1.2.1	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	Unique instance ID	ALWAYS	AUTO
Instance Creation Date	(0008,0012)	DA	Instance creation date	ALWAYS	AUTO
Instance Creation Time	(0008,0013)	TM	Instance creation time	ALWAYS	AUTO
Instance Number	(0020,0013)	IS		VNAP	AUTO

**Table 7.1-23  
PRESENTATION SERIES MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Modality	(0008,0060)	CS	PR	ALWAYS	AUTO
Presentation Creation Date	(0070,0082)	DA		ALWAYS	AUTO
Presentation	(0070,0083)	TM		ALWAYS	AUTO

# Planned

Creation Time					
Instance Number	(0020,0013)	IS		ALWAYS	AUTO
Content Label	(0070,0080)	CS	PR_ + image id	ALWAYS	AUTO
Content Description	(0070,0081)	LO		EMPTY	AUTO
Content Creator's Name	(0070,0084)	PN	PM	ALWAYS	AUTO

**Table 7.1-24  
PRESENTATION STATE IDENTIFICATION MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Presentation Creation Date	(0070,0082)	DA		ALWAYS	AUTO
Presentation Creation Time	(0070,0083)	TM		ALWAYS	AUTO
Instance Number	(0020,0013)	IS		ALWAYS	AUTO
Content Label	(0070,0080)	CS	PR_ + image id	ALWAYS	AUTO
Content Description	(0070,0081)	LO		EMPTY	AUTO
Content Creator's Name	(0070,0084)	PN	PM	ALWAYS	AUTO

**Table 7.1-25  
PRESENTATION STATE RELATIONSHIP MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Referenced Series Sequence	(0008,1115)	SQ		ALWAYS	AUTO
>Series Instance UID	(0020,000E)	UI		ALWAYS	AUTO
>Referenced Image Sequence	(0008,1140)	SQ		ALWAYS	AUTO
>>Referenced SOP Class UID	(0008,1150)	UI		ALWAYS	AUTO
>>Referenced SOP Instance UID	(0008,1155)	UI		ALWAYS	AUTO

**Table 7.1-26  
DISPLAYED AREA MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Displayed Area Selection Sequence	(0070,005A)	SQ		ALWAYS	AUTO

# Planned

>Referenced Image Sequence	(0008,1140)	SQ		ALWAYS	AUTO
>>Referenced SOP Class UID	(0008,1150)	UI		ALWAYS	AUTO
>>Referenced SOP Instance UID	(0008,1155)	UI		ALWAYS	AUTO
>Displayed Area Top Left Hand Corner	(0070,0052)	SL	0\0	ALWAYS	AUTO
>Displayed Area Bottom Right Hand Corner	(0070,0053)	SL	Image width and height	ALWAYS	AUTO
>Presentation Size Mode	(0070,0100)	CS	SCALE_TO_FIT	ALWAYS	AUTO
>Presentation Pixel Spacing	(0070,0101)	DS	The images pixel size	ALWAYS	AUTO

**Table 7.1-27  
GRAPHIC ANNOTATION MODULE OF CREATED SOP INSTANCES**

<b>Attribute Name</b>	<b>Tag</b>	<b>VR</b>	<b>Value</b>	<b>Presence of Value</b>	<b>Source</b>
Graphic Annotation Sequence	(0070,0001)	SQ		ALWAYS	AUTO
>Graphic Layer	(0070,0002)	CS	L + image id	ALWAYS	AUTO
>Text Object Sequence	(0070,0008)	SQ	Present if there are annotations	ANAP	AUTO
>>Bounding Box Annotation Units	(0070,0003)	CS	PIXEL	ALWAYS	AUTO
>>Unformatted Text Value	(0070,0006)	ST		ALWAYS	USER
>>Bounding Box Top Left Hand Corner	(0070,0010)	FL		ALWAYS	AUTO
>>Bounding Box Bottom Right Hand Corner	(0070,0011)	FL		ALWAYS	AUTO
>>Bounding Box Text Horizontal Justification	(0070,0012)	CS	LEFT	ALWAYS	AUTO
>Graphic Object Sequence	(0070,0009)	SQ	Present if there are graphic objects	ANAP	AUTO
>>Graphic	(0070,0005)	CS	PIXEL	ALWAYS	AUTO

# Planned

Annotation Units					
>>Graphic Dimensions	(0070,0020)	US	2	ALWAYS	AUTO
>>Number of Graphic Points	(0070,0021)	US		ALWAYS	AUTO
>>Graphic Data	(0070,0022)	FL	Graphic points data	ALWAYS	AUTO
>>Graphic Type	(0070,0023)	CS		ALWAYS	AUTO
>>Graphic Filled	(0070,0024)	CS	Y or N	ALWAYS	AUTO

**Table 7.1-28  
GRAPHIC LAYER SEQUENCE MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Graphic Layer Sequence	(0070,0060)	SQ		ALWAYS	AUTO
>Graphic Layer	(0070,0002)	CS	L + image id	ALWAYS	AUTO
Graphic Layer Order	(0070,0062)	IS	1	ALWAYS	AUTO
>Graphic Layer Recommended Display Grayscale Value	(0070,0066)	US	0xFFFF	ALWAYS	AUTO
>Graphic Layer Recommended Display RGB Value	(0070,0067)	US	0xFFFF\0xFFFF\0xFFFF (This attribute is retired and will be replaced by "Graphic Layer Recommended Display CIELab Value" in future software version)	ALWAYS	AUTO

**Table 7.1-29  
MODALITY LUT MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Rescale Intercept	(0028,1052)	DS	0	ALWAYS	AUTO
Rescale Slope	(0028,1053)	DS	1	ALWAYS	AUTO
Rescale Type	(0028,1054)	LO	US	ALWAYS	AUTO

**Table 7.1-30  
SOFTCOPY VOI LUT MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Softcopy VOI LUT	(0028,3110)	SQ		ALWAYS	AUTO

# Planned

Sequence					
Window Center	(0028,1050)	DS		ALWAYS	AUTO
Window Width	(0028,1051)	DS		ALWAYS	AUTO

**Table 7.1-31  
SOFTCOPY PRESENTATION LUT MODULE OF CREATED SOP INSTANCES**

Attribute Name	Tag	VR	Value	Presence of Value	Source
Presentation LUT Shape	(2050,0020)	CS		ALWAYS	AUTO

## 7.1.2 Usage of Attributes form received IOD's

The attributes Nuance Manager 3 uses from received IOD's are detailed in AE presentations above.

## 7.2 DATA DICTIONARY OF PRIVATE ATTRIBUTES

The Private Attributes added to created SOP Instances are listed in the Table Below. Nuance Manager 3 reserves a block of private attributes in the group 0009. Further details on usage of these private attributes are contained in Section 8.1.

**Table7.2-1  
DATA DICTIONARY OF PRIVATE ATTRIBUTES**

Tag	Attribute Name	VR	VM
(0101,0010)	Private Element Creator	LO	1
(0101,1005)	Image Enhancement Parameter File	LO	1
(0101,1006)	Original Sigmoid Ratio	IS	1

## 7.3 STANDARD EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES

Nuance Manager 3 does not implement any Specialized or Private SOP Classes.

## 7.4 PRIVATE TRANSFER SYNTAXES

Nuance Manager 3 does not implement any private Transfer Syntaxes.