



EN

Plannedmed Verity

IHE Integration Statement

Product name: Verity Manager

Version: 1.2

Date: June 2012

IHE Integration Statement			
Vendor	Product Name	Version	Date
Planmed	Verity Manager	1.2	June 2012
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Modality Group Case	
		Assisted Acquisition Protocol Setting	
Patient Information Reconciliation	Acquisition Modality	none	
Internet address for vendor's IHE information: http://www.planmed.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	not applicable		
DICOM	http://www.planmed.com/dicom		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	